

Product Datasheet

**NGF Receptor Antibody / p75NTR / CD271
(orb534446)**

Description

CD271/NGFR pays a role in the

to the cell surface in adipoc...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

A portion of amino acids 281-421 from the human protein was used as the immunogen for the recombinant NGFR antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the recombinant NGFR antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The stated application concentrations are suggested starting points. Titration of the recombinant NGFR antibody may be required due to differences in protocols and secondary/substrate sensitivity.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P08138

Hazard
Information

This recombinant NGFR antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml