

Product Datasheet

PLIN2 antibody (orb534564)

Description

Rabbit monoclonal antibody to PLIN2

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications
Immunogen

Recombinant fragment (around aa 249-376) of human Adipophilin (ADFP) protein (exact sequence is proprietary)

Target

PLIN2

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Rabbit / IgG

Clonality

Monoclonal

MW

48kDa

Uniprot ID
[Q99541](#)
Entrez
[123](#)
Dilution Range

WB: 1-2 µg/ml , IHC-P: 1-2 µg/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.