

Product Datasheet

Elastin antibody (orb534594)

Description	Rabbit monoclonal antibody to Elastin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human Elastin protein
Target	ELN
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	70kDa
Uniprot ID	P15502
Entrez	2006
Dilution Range	IHC-P: 1-2 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.