

Product Datasheet

Ferritin antibody (orb534629)

Description

Mouse monoclonal antibody to Ferritin

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant fragment (around aa 38-165) of human FTL protein (exact sequence is proprietary)
Target	FTL
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	19-25kDa
Uniprot ID	P02792
Entrez	2512
Dilution Range	WB: 0.5-2 µg/ml, IHC-P: 0.1-0.2 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.