

Product Datasheet

GFAP antibody (orb534640)

Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Bovine, Gallus, Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant full-length human GFAP protein
Target	GFAP
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	~50kDa
Uniprot ID	P14136
Entrez	2670
Dilution Range	WB: 0.5-1 µg/ml , IHC-P: 0.25-0.5 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.