

Product Datasheet

GP2 antibody (orb534647)

Description

Mouse monoclonal antibody to GP2

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications
Immunogen

Recombinant fragment of human GP2 protein (around aa 35-179) (exact sequence is proprietary)

Target

GP2

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG2c, kappa

Clonality

Monoclonal

MW

59kDa

Uniprot ID

P55259

Entrez

2813

Dilution Range

ELISA: 2-4 µg/ml, IHC-P: 1-2 µg/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.