


Product Datasheet

LRRK1 antibody (orb53467)

Description	LRRK1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Leucine-rich repeat serine/threonine-protein kinase 1 protein (1-68AA)
Target	LRRK1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q38SD2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2(lane 1), ...