

Product Datasheet

IL3RA antibody (orb534697)

Description	Rabbit monoclonal antibody to IL3RA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Recombinant fragment of human IL3RA protein (around aa 26-171) (exact sequence is proprietary)
Target	IL3RA
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	~70kDa
Uniprot ID	P26951
Entrez	3563
Dilution Range	IHC-P: 0.25-0.5 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.