

Product Datasheet

IVL antibody (orb534707)

Description	Rabbit monoclonal antibody to IVL
Species/Host	Rabbit
Reactivity	Canine, Human, Monkey, Porcine, Primate
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human full-length Involucrin protein (IVRN)
Target	IVL
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	66-170kDa
Uniprot ID	P07476
Entrez	3713
Dilution Range	FACS: 0.5-1 µg/million cells , IF: 1-2 µg/ml , IHC-P: 0.1-0.2 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.