

Product Datasheet

KRT6A antibody (orb534711)

Description	Mouse monoclonal antibody to KRT6A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	Recombinant full-length human Cytokeratin 6A (KRT6A)protein
Target	KRT6A
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	56kDa
Uniprot ID	P02538
Entrez	3853
Dilution Range	IHC-P: 1-2 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.