

Product Datasheet

KRT18 antibody (orb534724)

Description

Rabbit monoclonal antibody to KRT18

Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Feline, Ferret, Goat, Hamster, Human, Monkey, Mouse, Porcine, Rat, Sheep

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

Recombinant human full-length KRT18 protein

Target

KRT18

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Rabbit / IgG

Clonality

Monoclonal

MW

45kDa

Uniprot ID
P05783
Entrez
3875
Dilution Range

WB: 0.5-2 µg/ml, IHC-P: 0.5-1 µg/ml for 30 min at RT

Expiration Date

12 months from date of receipt.