

Product Datasheet

MKI67 antibody (orb534748)

Description

Mouse monoclonal antibody to MKI67

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant human Ki67 protein fragment (around aa 2293-2478) (exact sequence is proprietary)
Target	MKI67
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b, kappa
Clonality	Monoclonal
MW	345kDa and 395kDa
Uniprot ID	P46013
Entrez	4288
Dilution Range	FACS: 1-2 µg/million cells, IF: 1-2 µg/ml , IHC-P: 1-2 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.