

Product Datasheet

MKI67 antibody (orb534751)

Description

Mouse monoclonal antibody to MKI67

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P

Applications
Immunogen

Recombinant fragment of human Ki67 protein (around aa 2293-2478) (exact sequence is proprietary)

Target

MKI67

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG2b, kappa

Clonality

Monoclonal

MW

345kDa and 395kDa

Uniprot ID

P46013

Entrez

4288

Dilution Range

FACS: 1-2 µg/million cells, IF: 1-2 µg/ml, IHC-P: 1-2 µg/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.