

Product Datasheet

MYH11 antibody (orb534773)

Description	Rabbit monoclonal antibody to MYH11
Species/Host	Rabbit
Reactivity	Bovine, Canine, Feline, Gallus, Guinea pig, Human, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human MYH11 protein
Target	MYH11
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	205kDa (MHC-1) and 200kDa (MHC-2)
Uniprot ID	P35749
Entrez	4629
Dilution Range	IHC-P: 0.5-1 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.