

Product Datasheet

NGFR antibody (orb534782)

Description	Mouse monoclonal antibody to NGFR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant fragment of human p75 NGFR protein (around aa 281-421) (exact sequence is proprietary)
Target	NGFR
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b, kappa
Clonality	Monoclonal
MW	75kDa
Uniprot ID	P08138
Entrez	4804
Dilution Range	IHC-P: 1-2 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.