

Product Datasheet

TIMP1 antibody (orb534916)

Description	Rabbit monoclonal antibody to TIMP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant human full-length TIMP1 protein
Target	TIMP1
Preservatives	200ug/ml of Ab Purified by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	23kDa
Uniprot ID	P01033
Entrez	7076
Dilution Range	FACS: 0.5-1 µg/million cells , WB: 0.5-1.0 µg/ml , IHC-P: 0.5-1.0 µg/ml for 30 minutes at RT, IF: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.