

Product Datasheet

TIMP2 antibody (orb534918)

Description	Rabbit monoclonal antibody to TIMP2
Species/Host	Rabbit
Reactivity	Bovine, Guinea pig, Human, Mouse, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A synthetic peptide from the N-terminus of human TIMP2 protein (aa56-66) (DSGNDIYGNPI)
Target	TIMP2
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	21kDa
Uniprot ID	P16035
Entrez	7077
Dilution Range	IHC-P: 1-2 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.