

Product Datasheet

p53 antibody (orb534930)

Description	Rabbit monoclonal antibody to p53
Species/Host	Rabbit
Reactivity	Canine, Gallus, Hamster, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human TP53 protein
Target	TP53
Preservatives	200ug/ml of Ab Purified by Protein A Column. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	53kDa
Uniprot ID	P04637
Entrez	7157
Dilution Range	WB: 0.5-1 µg/ml , IHC-P: 0.5-1 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.