

Product Datasheet

Ezrin antibody (orb534949)

Description	Mouse monoclonal antibody to Ezrin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human full-length Ezrin protein
Target	Ezrin
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b
Clonality	Monoclonal
MW	87kDa
Uniprot ID	P15311
Entrez	7430
Dilution Range	WB: 0.5-1.0 µg/ml , IHC-P: 0.5-1 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.