

Product Datasheet

CCND2 antibody (orb534978)

Description	Mouse monoclonal antibody to CCND2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant full-length human Cyclin D2 (CCND2) protein
Target	CCND2
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG, kappa
Clonality	Monoclonal
MW	34kDa
Uniprot ID	P30279
Entrez	894
Dilution Range	IHC-P: 1-2 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.