

Product Datasheet

CDC2 antibody (orb535045)

Description	Mouse monoclonal antibody to CDC2
Species/Host	Mouse
Reactivity	Gallus, Guinea pig, Human, Mouse, Rat, Xenopus
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, IP, WB
Immunogen	C-Terminal 2/3rds of Xenopus cdc2 expressed in E.coli
Target	CDC2
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	34kDa
Uniprot ID	P06493
Entrez	983
Dilution Range	IF: 1-2 µg/ml , IP: 2-4 µg/mg protein, WB: 0.5-1.0 µg/ml, IHC-P: 2-4 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.