

Product Datasheet

CDH1 antibody (orb535051)

Description	Rabbit monoclonal antibody to CDH1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant full-length human E-Cadherin protein
Target	CDH1
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	120-80kDa (Mature); 135kDa (Precursor)
Uniprot ID	P12830
Entrez	999
Dilution Range	WB: 0.5-1 µg/ml , IHC-P: 1-2 µg/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.