

Product Datasheet

CYP20A1 antibody (Biotin) (orb53506)

Description	CYP20A1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cytochrome P450 20A1 protein (91-390AA)
Target	CYP20A1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6UW02
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.