

Product Datasheet

KRT77 & KRT76 antibody (orb535060)

Description	Rabbit monoclonal antibody to KRT77 & KRT76
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human KRT76 and KRT77 proteins
Target	KRT77 & KRT76
Preservatives	200ug/ml of recombinant MAb Purified by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	40-67kDa
Uniprot ID	Q01546 , Q7Z794
Entrez	374454 , 51350
Dilution Range	IHC-P: 0.25-0.5 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.