

Product Datasheet

Cytokeratin antibody (orb535062)

Description

Rabbit monoclonal antibody to Cytokeratin

Species/Host

Rabbit

Reactivity

Bovine, Frog, Goat, Guinea pig, Human, Monkey, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

Keratin-enriched preparation from cultured human epithelial cells

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C.
Antibody without azide - store at -20 to -80°C.
Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Rabbit / IgG

Clonality

Monoclonal

MW

Multiple

Uniprot ID

P04259, P19013, P02538, P48668, P13647

Entrez

3860, 3875, 3851, 3852, 3856, 286887, 3858, 3853, 3854

Dilution Range

FACS: 0.5-1 µg/million cells , IF: 0.5-1.0 µg/ml
WB: 0.5-1.0 µg/ml , IHC-P: 0.5-1.0 µg/ml for 30 min at RT

Expiration Date

12 months from date of receipt.