

## Product Datasheet

### Bovine COL1A1 protein (orb535090)

<b>Description</b>	Bovine COL1A1 protein
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	None
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	Store vial at 4° C prior to opening. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration of 6 months after opening if stored, handled and kept under sterile conditions.
<b>Note</b>	For research use only
<b>Application notes</b>	Ultrapure Bovine Collagen Type I solution are suitable for plate coating. Enough to coat 10 (ten) T75 plates.
<b>Purity</b>	Ultrapure Bovine Collagen Type I Solution is purified from bovine placenta and shown to be 99% pure by gel electrophoresis and immunoassays.
<b>Expiration Date</b>	Please enquire.