

Product Datasheet

Bovine COL1A1 protein (orb535091)

Description	Recombinant of Bovine COL1A1 protein
Conjugation	Unconjugated
Target	COL1A1
Preservatives	Preservative: None. Stabilizer: None. 0.013M HCl
Form/Appearance	Liquid (sterile filtered)
Concentration	2.0mg/ml
Storage	Store vial at 4° C prior to opening. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration of 6 months after opening if stored, handled and kept under sterile conditions.
Note	For research use only
Application notes	Ultrapure Bovine Collagen Type I solution has been tested in cell culture and can be utilized to mimic in vivo environment for 2D and 3D cell culture, tissue engineering research, and biochemistry.
Purity	Ultrapure Bovine Collagen Type I Solution is purified from bovine placenta and shown to be greater than 99% pure by gel electrophoresis and immunoassays.
Expiration Date	6 months from date of receipt.