

Product Datasheet

12 Epitope GST Tag Protein Marker Lysate (orb535115)

Description

12 Epitope GST Tag Protein Marker

Conjugation

Unconjugated

Tested

SDS-PAGE, WB

Applications
Preservatives

None. 1X SDS-PAGE Sample Buffer
(62.5 mM Tris HCl, 2% SDS, 10%
Glycerol and 0.005% bromophenol
blue, pH 6.8)

Form/Appearance

Liquid

Concentration

5 mg/ ml

Storage

Store vial at -20° C prior to opening.
Aliquot contents and freeze at -20° C
or below for extended storage. Avoid
cycles of freezing and thawing.

Note

For research use only

Application notes

This 12 Epitope Tag Protein Marker Lysate has been tested in SDS-PAGE and western blot and is designed for use as a positive control for anti-epitope tag antibodies. This E.coli lysate contains a recombinant protein comprised of the following commonly used epitope tags: GST--T7--HSV--cMyc--VSV--Glu-Glu--V5--E-tag--Flag--S tag--HA--6XHis. A 48 kDa recombinant protein is visualized after SDS-polyacrylamide gel electrophoresis and western blotting using an antibody recognizing any of the expressed epitope tags listed above. Thaw the lysate to room temperature. Gently invert vial several times to ensure that the solution is homogeneous and that any precipitated material is re-dissolved. A loading volume of 0.25 to 2 µL is recommended for western blotting.

Purity

Loading control lysate made in E. Coli for Western blotting expresses 12 different epitope tags.

Hazard

Non-Toxic

Information
Dilution Range

WB: User Optimized

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558