

Product Datasheet

CBP antibody (orb535293)

Description	Rabbit polyclonal antibody against CREBBP conjugated to HRP
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the mid-protein of Human CBP (aa. 2100-2200)
Target	CBP
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	#REF!
Clonality	Polyclonal
MW	265 kDa
Uniprot ID	Q92793
NCBI	NP_001073315.1
Entrez	1387
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.