

## Product Datasheet

### CBP antibody (orb535295)

|                            |                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody against CREBBP conjugated to PerCP                              |
| <b>Species/Host</b>        | Rabbit                                                                                     |
| <b>Reactivity</b>          | Human                                                                                      |
| <b>Conjugation</b>         | PerCP                                                                                      |
| <b>Tested Applications</b> | IHC, WB                                                                                    |
| <b>Immunogen</b>           | Synthetic peptide from the mid-protein of Human CBP (aa. 2100-2200)                        |
| <b>Target</b>              | CBP                                                                                        |
| <b>Preservatives</b>       | PBS pH 7.4, 50% glycerol, 0.09% sodium azide<br>*Storage buffer may change when conjugated |
| <b>Concentration</b>       | 1 mg/ml                                                                                    |
| <b>Storage</b>             | Conjugated antibodies should be stored at 4°C                                              |
| <b>Note</b>                | For research use only                                                                      |
| <b>Application notes</b>   | #REF!                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                 |
| <b>MW</b>                  | 265 kDa                                                                                    |
| <b>Uniprot ID</b>          | <a href="#">Q92793</a>                                                                     |
| <b>NCBI</b>                | <a href="#">NP_001073315.1</a>                                                             |
| <b>Entrez</b>              | <a href="#">1387</a>                                                                       |
| <b>Dilution Range</b>      | WB (1:1000); IHC (1:50)                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                            |