

Product Datasheet

ROR1 antibody (orb535471)

Description	Rabbit polyclonal antibody against ROR1 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the mid-protein of Human ROR1 (aa. 200-300)
Target	ROR1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	#REF!
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q01973
NCBI	NP_001077061.1
Entrez	4919
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.