

Product Datasheet

KRAS antibody (orb535740)

Description	Rabbit polyclonal antibody against KRAS conjugated to APC
Species/Host	Rabbit
Reactivity	Rat
Conjugation	APC
Tested Applications	WB
Immunogen	Synthetic peptide from the C-terminal of Human KRAS (aa. 100-189)
Target	KRAS
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	#REF!
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	P01116
NCBI	NP_004976.2
Entrez	3845
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.