

Product Datasheet

MDM2 antibody [Out of stock] (orb535749)

Description	Rabbit polyclonal antibody against MDM2 conjugated to ATTO-390
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-390
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the mid-protein of Human MDM2 (aa. 100-200)
Target	MDM2
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	#REF!
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q00987
NCBI	NP_001138811.1
Entrez	4193
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.