

Product Datasheet

PKC gamma antibody [Out of stock] (orb535822)

Description

Rabbit polyclonal antibody against

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-488
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the mid-protein of Human PRKCG (aa. 150-250)
Target	PRKCG
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-782 was sufficient for detection of PRKCG in 15 µg of mouse Brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	62 kDa, 40 kDa
Uniprot ID	P05129
NCBI	NP_001303258.1
Entrez	5582
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.