

Product Datasheet

SFRP2 antibody (orb535903)

Description	Rabbit polyclonal antibody against SFRP2 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the C-terminal of Human Secreted frizzled-related protein 2 (aa. 195-295)
Target	SFRP2
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-786 was sufficient for detection of SFRP2 in 15 µg of mouse Brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	Q96HF1
NCBI	NP_003004.1
Entrez	6423
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.