

Product Datasheet

DKK1 antibody [Out of stock] (orb535912)

Description

Rabbit polyclonal antibody against DKK1

Species/Host	Rabbit
Reactivity	Human
Conjugation	ATTO-488
Tested Applications	WB
Immunogen	Synthetic peptide from the mid-protein of Human Dickkopf-related protein 1 (aa. 120-220)
Target	DKK1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-787 was sufficient for detection of DKK1 in 10 µg of human 293T cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	O94907
NCBI	NP_036374.1
Entrez	22943
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.