

Product Datasheet

STK4 antibody (orb536038)

Description	Rabbit polyclonal antibody against STK4 conjugated to ATTO-390
Species/Host	Rabbit
Reactivity	Human
Conjugation	ATTO-488
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the C-terminal of Human Serine/threonine-protein kinase 4 (aa. 320-420)
Target	STK4
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-794 was sufficient for detection of STK4 in 10 µg of human HeLa cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	Q13043
NCBI	NP_006273.1
Entrez	6789
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.