

Product Datasheet

RICTOR antibody (orb536106)

Description

Rabbit polyclonal antibody against RICTOR

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	RPE
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the mid-protein of Human Rapamycin-insensitive companion of mTOR (aa. 1060-1160)
Target	RICTOR
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-797 was sufficient for detection of RICTOR in 15 µg of mouse Brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	192 kDa
Uniprot ID	Q6R327
NCBI	NP_001272368.1
Entrez	253260
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.