

Product Datasheet

**Alpha Synuclein antibody [Out of stock]
(orb536146)**

Description

Rabbit polyclonal antibody against

390

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-488

Tested

ELISA, IHC, WB

Applications
Immunogen

Recombinant Human alpha Synuclein Fibril

Target

Alpha Synuclein

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-800 was sufficient for detection of Alpha Synuclein in 15 µg of mouse Brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

14.46kDa

Uniprot ID
P37840
NCBI
NP_000336.1
Entrez
6622
Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.