

Product Datasheet

Alpha Synuclein antibody (orb536155)

Description

Rabbit polyclonal antibody against alpha

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human alpha Synuclein Fibril
Target	Alpha Synuclein
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-800 was sufficient for detection of Alpha Synuclein in 15 µg of mouse Brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	14.46kDa
Uniprot ID	P37840
NCBI	NP_000336.1
Entrez	6622
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.