

Product Datasheet

TSR2 antibody (HRP) (orb53616)

Description	TSR2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Pre-rRNA-processing protein TSR2 homolog protein (1-191AA)
Target	TSR2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969E8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.