

Product Datasheet

Tau antibody [Out of stock] (orb536166)

Description	Rabbit polyclonal antibody against Tau conjugated to ATTO-390
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	ATTO-594
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Tau 2NR4 P301S Fibril
Target	Tau
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-801 was sufficient for detection of Tau in 10 µg of human SH-SY5Y cell lysate cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	45.8kDa
Uniprot ID	P10636
Entrez	4137
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.