

Product Datasheet

MNT antibody (orb536204)

Description	Goat polyclonal antibody to MNT
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ROX / MNT
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70-75kDa band observed in Mouse Brain lysates (calculated MW of 63.3kDa according to NP_034943.3). This molecular weight is routinely observed by other sources. Recommended concentration: 0.3-1µg/ml. Primary incubation 1 hour at room temperature.
Protein Sequence	DQASTSTASEGEDNI
Clonality	Polyclonal
MW	62.3
Entrez	4335
Dilution Range	ELISA: 1:4000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.