

Product Datasheet

CTSB antibody (orb536205)

Description	Goat polyclonal antibody to CTSB
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	Cathepsin B (aa194-206)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the vesicles of HeLa cells. Recommended concentration: 10µg/ml.
Protein Sequence	SRPPCTGEGDTPK
Clonality	Polyclonal
MW	37.8; 23.7
Entrez	1508
Dilution Range	IF: 10 µg/ml, P-ELISA: 1:2000
Expiration Date	12 months from date of receipt.