


Product Datasheet

CXXC1 antibody (orb536270)

Description	Rabbit polyclonal antibody to CXXC1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human CXXC1(NP_001095124.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80kDa
Uniprot ID	Q9P0U4
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A-549 cells...