


Product Datasheet

DCP1A antibody (orb536271)

Description	Rabbit polyclonal antibody to DCP1A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human DCP1A (NP_060873.4).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	70KDa
Uniprot ID	Q9NPI6
Entrez	55802
Dilution Range	IF: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...