


Product Datasheet

FABP6 antibody (orb536273)

Description	Rabbit polyclonal antibody to FABP6
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human FABP6(NP_001436.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa/19kDa
Uniprot ID	P51161
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...