

Product Datasheet

IP6K2 antibody (orb536278)

Description	Rabbit polyclonal antibody to IP6K2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human IP6K2(NP_057375.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8kDa/9kDa/11kDa/20kDa/21kDa/49kDa
Uniprot ID	Q9UHH9
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.