


Product Datasheet

LHX6 antibody (orb536284)

Description	Rabbit polyclonal antibody to LHX6
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human LHX6(NP_055183.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19kDa/38kDa/40kDa/41kDa/43kDa
Uniprot ID	Q9UPM6
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...