


Product Datasheet

POMGNT2 antibody (orb536294)

Description	Rabbit polyclonal antibody to POMGNT2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human POMGNT2(NP_116195.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66kDa
Uniprot ID	Q8NAT1
Dilution Range	IF: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...